


Anti-Usp9X Polyclonal Antibody For Wb, Elisa - Q93008

Item Number: CM0022865



Introduction

USP9X Rabbit Polyclonal Antibody for Western Blot and ELISA applications, validated in Human, Mouse, Rat. High-quality IVD raw material for studying deubiquitinase activity, TGF- β signaling, and circadian rhythms.

[Learn More](#)

Parameter	Value
Product Name	USP9X Rabbit pAb
Remarks/Alias	FAF; FAM; hFAM; DFFRX; MRX99; XLID99; MRXS99F; USP9X
Species	Human
Gene ID (Human)	8239
Gene ID	8239
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human USP9X (NP_001034679.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93008
Protein Weight	290kDa
Shipping	Ice bag