

## Anti-Usp9X/Usp9Y Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00507

Item Number: CM0023038

### Introduction

Rabbit polyclonal antibody targeting human USP9X/USP9Y. Validated for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Suitable for deubiquitination and male germ cell development studies. Molecular weight approximately 291kDa.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | USP9X/USP9Y Rabbit pAb                                                                                           |
| Remarks/Alias    | DFFRY; SPGFY2; USP9X/USP9Y                                                                                       |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 8287                                                                                                             |
| Gene ID          | 8287                                                                                                             |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 2456-2555 of human USP9X/USP9Y (NP_004645.2). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | O00507                                                                                                           |
| Protein Weight   | 291kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                          |