


Anti-Utp18 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y5J1

Item Number: CM0020818



Introduction

High-quality rabbit polyclonal antibody targeting human UTP18 (WDR50), validated for Western blot and ELISA. Ideal for investigating ribosome biogenesis, SSU processome, and pre-18S rRNA processing.

[Learn More](#)

Parameter	Value
Product Name	UTP18 Rabbit pAb
Remarks/Alias	WDR50; CGI-48; UTP18
Species	Human
Gene ID (Human)	51096
Gene ID	51096
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 70-260 of human UTP18 (NP_057085.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5J1
Protein Weight	62kDa
Shipping	Ice bag