


Anti-Utp6 Polyclonal Antibody For Wb, Elisa - Q9Nyh9

Item Number: CM0019179



Introduction

Rabbit polyclonal antibody against human UTP6 (HCA66), validated for Western Blot and ELISA. Cross-reacts with mouse and rat. Ideal for ribosome biogenesis and nucleolar processing studies.

[Learn More](#)

Parameter	Value
Product Name	UTP6 Rabbit pAb
Remarks/Alias	HCA66; C17orf40; UTP6
Species	Human
GeneID (Human)	55813
GeneID	55813
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYH9
Protein Weight	70kDa
Shipping	Ice bag