

Anti-Uvrag Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P2Y5

Item Number: CM0027727



Introduction

Rabbit polyclonal antibody against human UVRAG (p63, VPS38). Validated for Western blot, IHC-P, and ELISA applications. Detects endogenous UVRAG in human, mouse, and rat samples. Suitable for autophagy, membrane trafficking, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	UVRAG Rabbit pAb
Remarks/Alias	p63; DHTX; VPS38; UVRAG
Species	Human
GeneID (Human)	7405
GeneID	7405
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 501-699 of human UVRAG (NP_003360.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2Y5
Protein Weight	78kDa
Shipping	Ice bag