

Anti-B-Amyloid Monoclonal Antibody For If-P, Ihc-P, Db, Elisa - P05067

Item Number: CM0010487



Introduction

Rabbit monoclonal antibody to human β -Amyloid (P05067). Validated for IF-P, IHC-P, DB, ELISA; cross-reacts with mouse. Suitable for Alzheimer's disease research and amyloid plaque analysis.

[Learn More](#)

Parameter	Value
Product Name	β -Amyloid Rabbit mAb
Remarks/Alias	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP
Species	Human
GeneID (Human)	351
GeneID	351
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, DB, ELISA
Cross Reactivity	Mouse
SWISS	P05067
Protein Weight	87kDa
Shipping	Ice bag