

Anti-B-Amyloid(1-40) Rabbit Monoclonal Antibody For If-P, Ihc-P, Db, Elisa - P05067

Item Number: CM0006297



Introduction

High-affinity rabbit monoclonal antibody targeting β -Amyloid(1-40), a key Alzheimer's disease peptide. Validated for IF-P, IHC-P, DB, ELISA. Cross-reacts with mouse. Ideal for amyloid and APP processing studies.

[Learn More](#)

Parameter	Value
Product Name	β -Amyloid(1-40) Rabbit mAb
Remarks/Alias	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP; β -Amyloid(1-40)
Species	Human
GeneID (Human)	351
GeneID	351
Immunogen	A synthetic peptide corresponding to amino acids 671-770 of human β -Amyloid(1-40) (NP_000475.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, DB, ELISA
Cross Reactivity	Mouse
SWISS	P05067
Protein Weight	4kDa
Shipping	Ice bag