

Anti-B-Arrestin1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P49407

Item Number: CM0013360



Introduction

Rabbit polyclonal antibody targeting β -arrestin1 (ARRB1/ARB1/ARR1), validated in WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Suitable for GPCR desensitization and internalization studies.

[Learn More](#)

Parameter	Value
Product Name	β -arrestin1 Rabbit pAb
Remarks/Alias	ARB1; ARR1; β -arrestin1
Species	Human
Gene ID (Human)	408
Gene ID	408
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 169-418 of human β -arrestin1 (NP_004032.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49407
Protein Weight	47kDa
Shipping	Ice bag