

Anti-B-Catenin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35222

Item Number: CM0018015



Introduction

β -Catenin (CTNNB1) rabbit polyclonal antibody validated for Western blot, IHC, IF, and ELISA. Reacts with human, mouse, and rat samples. Ideal for Wnt signaling and cell adhesion research.

[Learn More](#)

Parameter	Specification
Product Name	β -Catenin Rabbit pAb
Remarks/Alias	EVR7; CTNNB; MRD19; NEDSDV; armadillo; β -Catenin
Species	Human
GeneID (human)	1499
GeneID	1499
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human beta Catenin (NP_001895.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35222
Protein Weight	85kDa
Shipping	Ice bag