

## Anti-B-Catenin Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P35222

Item Number: CM0015102



### Introduction

Rabbit polyclonal antibody targeting human  $\beta$ -Catenin (CTNNB1), suitable for WB, IHC-P, IP, and ELISA applications. Detects endogenous  $\beta$ -Catenin from human, mouse, and rat samples. Ideal for Wnt signaling and cell adhesion research.

[Learn More](#)

| Parameter        | Value                                                                                                                                             |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | $\beta$ -Catenin Rabbit pAb                                                                                                                       |
| Remarks/Alias    | EVR7; CTNNB; MRD19; NEDSDV; armadillo; $\beta$ -Catenin                                                                                           |
| Species          | Human                                                                                                                                             |
| Gene ID (Human)  | 1499                                                                                                                                              |
| Gene ID          | 1499                                                                                                                                              |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 667-781 of human $\beta$ -Catenin (NP_001895.1) |
| Source           | Rabbit                                                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                                                             |
| Application      | WB, IHC-P, IP, ELISA                                                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                                 |
| SWISS            | P35222                                                                                                                                            |
| Protein Weight   | 85kDa                                                                                                                                             |
| Shipping         | Ice bag                                                                                                                                           |