

Anti-B-Trcp/Btrc Rabbit Mab For Wb, Ihc-P, Elisa - Q9Y297

Item Number: CM0006350



Introduction

Recombinant monoclonal antibody targeting human β -TrCP/BTRC, validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat BTRC. Ideal for ubiquitination and Wnt signaling studies.

[Learn More](#)

Parameter	Value
Product Name	β -TrCP/BTRC Rabbit mAb
Remarks/Alias	FWD1; FBW1A; FBXW1; bTrCP; FBXW1A; bTrCP1; betaTrCP; BETA-TRCP; β -TrCP/BTRC
Species	Human
GeneID (human)	8945
GeneID	8945
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human β -TrCP/BTRC (NP_378663.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y297
Protein Weight	69kDa
Shipping	Ice bag