

Anti-B-Trcp/Btrc Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y297

Item Number: CM0022862



Introduction

Rabbit polyclonal antibody against human β -TrCP/BTRC (FWD1, FBXW1) for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for studying ubiquitin-mediated proteolysis, NF- κ B and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	β -TrCP/BTRC Rabbit pAb
Remarks/Alias	FWD1; FBW1A; FBXW1; bTrCP; FBXW1A; bTrCP1; betaTrCP; BETA-TRCP; β -TrCP/BTRC
Species	Human
GeneID	8945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human β -TrCP/BTRC (NP_378663.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y297
Protein Weight	69kDa
Shipping	Ice bag