

Anti-Vamp1 Monoclonal Antibody For Wb - P23763

Item Number: CM0004954

Introduction

Rabbit monoclonal antibody against human VAMP1/synaptobrevin-1, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. For neuroscience and vesicle trafficking research.

[Learn More](#)

Parameter	Value
Product Name	VAMP1 Rabbit mAb
Remarks/Alias	SAX1; SYB1; CMS25; SPAX1; VAMP-1; VAMP1
Species	Human
Gene ID	6843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human VAMP1 (P23763)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23763
Protein Weight	13kDa
Shipping	Ice bag