

Anti-Vamp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P23763

Item Number: CM0025503

Introduction

High-quality rabbit polyclonal antibody targeting VAMP1 (synaptobrevin-1) for Western blot, immunohistochemistry, and ELISA applications. Reacts with human, mouse, and rat VAMP1 (predicted ~13 kDa). Supports vesicle trafficking and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	VAMP1 Rabbit pAb
Remarks/Alias	SAX1; SYB1; CMS25; SPAX1; VAMP-1; VAMP1
Species	Human
Gene ID (Human)	6843
Gene ID	6843
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-96 of human VAMP1 (NP_954740.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23763
Protein Weight	13kDa
Shipping	Ice bag