


Anti-Vamp2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P63027

Item Number: CM0008946



Introduction

High-quality Anti-VAMP2 rabbit monoclonal antibody for WB, IHC-P, IF-P, ELISA applications. Reacts with human, mouse, rat VAMP2. Suitable for research on synaptic vesicle exocytosis and SNARE-mediated membrane fusion.

[Learn More](#)

Parameter	Value
Product Name	VAMP2 Rabbit mAb
Remarks/Alias	SYB2; VAMP-2; NEDHAHM; VAMP2
Species	Human
Gene ID	6844
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VAMP2 (P63027)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63027
Protein Weight	13kDa
Shipping	Ice bag