

Anti-Vamp3/Cellubrevin Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q15836

Item Number: CM0009332

Introduction

Rabbit monoclonal antibody against human VAMP3/cellubrevin, validated for WB, IF/ICC, and ELISA. Ideal for studying endosomal trafficking and SNARE-mediated vesicular transport in research settings.

[Learn More](#)

Parameter	Value
Product Name	VAMP3/Cellubrevin Rabbit mAb
Remarks/Alias	CEB; VAMP3/Cellubrevin
Species	Human
GenelD (human)	9341
GenelD	9341
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human VAMP3/Cellubrevin (Q15836).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q15836
Protein Weight	11kDa
Shipping	Ice bag