

Anti-Vamp4 Rabbit Polyclonal Antibody For Wb And Elisa - 075379

Item Number: CM0024962

Introduction

Rabbit anti-VAMP4 polyclonal antibody for Western blot (WB) and ELISA applications. Detects human, mouse, and rat VAMP4, a v-SNARE regulating vesicle fusion and membrane trafficking. Ideal for SNARE complex and muscle regeneration research.

[Learn More](#)

Parameter	Value
Product Name	VAMP4 Rabbit pAb
Remarks/Alias	VAMP-4; VAMP24; VAMP4
Species	Human
Gene ID (Human)	8674
Gene ID	8674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human VAMP4 (NP_003753.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75379
Protein Weight	16kDa
Shipping	Ice bag