


Anti-Vamp7 Polyclonal Antibody For Wb, If/Icc, Elisa - P51809

Item Number: CM0022914



Introduction

High-quality rabbit polyclonal antibody against human VAMP7 (SYBL1/TI-VAMP), validated for Western blot, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Suitable for lysosomal trafficking and exocytosis research.

[Learn More](#)

Parameter	Value
Product Name	VAMP7 Rabbit pAb
Remarks/Alias	SYBL1; TIVAMP; VAMP-7; TI-VAMP; VAMP7
Species	Human
Gene ID (Human)	6845
Gene ID	6845
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-70 of human VAMP7 (NP_005629.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51809
Protein Weight	25kDa
Shipping	Ice bag