

Anti-Vamp8 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bv40

Item Number: CM0014650



Introduction

The VAMP8 Rabbit Polyclonal Antibody (WB, ELISA) targets human VAMP8 (also known as Endobrevin/EDB) with cross-reactivity to mouse and rat. It is ideal for studying SNARE-mediated membrane fusion, autophagy, and immune cell exocytosis.

[Learn More](#)

Parameter	Value
Product Name	VAMP8 Rabbit pAb
Remarks/Alias	EDB; VAMP-8; VAMP8
Species	Human
Gene ID (Human)	8673
Gene ID	8673
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-75 of human VAMP8 (NP_003752.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV40
Protein Weight	11kDa
Shipping	Ice bag