


Anti-Vapa Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P0L0

Item Number: CM0016356



Introduction

Rabbit polyclonal antibody against human VAPA (VAP-A), validated for Western blot, IHC-P, and ELISA. Detects endogenous VAPA in human and mouse samples. Ideal for ER contact site and vesicle trafficking research.

[Learn More](#)

Parameter	Specification
Product Name	VAPA Rabbit pAb
Remarks/Alias	VAP-A; VAP33; VAMP-A; VAP-33; hVAP-33; VAPA
Species	Human
GeneID (human)	9218
GeneID	9218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 137-227 of human VAPA (NP_919415.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P0L0
Protein Weight	28kDa
Shipping	Ice bag