

Anti-Vasp Rabbit Monoclonal Antibody For Wb, Elisa - P50552

Item Number: CM0011846



Introduction

Rabbit monoclonal antibody against human VASP, validated for WB and ELISA. Ideal for studying actin dynamics, platelet activation, and cell migration. Detects a 40 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	VASP Rabbit mAb
Remarks/Alias	VASP; vasodilator-stimulated phosphoprotein
Species	Human
Gene ID (Human)	7408
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human VASP (P50552)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P50552
Protein Weight	40kDa
Shipping	Ice bag