

Anti-Vasp Rabbit Polyclonal Antibody For Wb, Elisa - P50552

Item Number: CM0016310



Introduction

VASP Rabbit pAb, a polyclonal antibody validated for WB and ELISA, specifically detects human vasodilator-stimulated phosphoprotein. Ideal for cytoskeleton dynamics and platelet activation research.

[Learn More](#)

Parameter	Specification
Product Name	VASP Rabbit pAb
Remarks/Alias	VASP
Species	Human
GeneID (Human)	7408
GeneID	7408
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 150-250 of human VASP (NP_003361.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P50552
Protein Weight	40kDa
Shipping	Ice bag