

Anti-Vasp Rabbit Polyclonal Antibody - P50552

Item Number: CM0015398



Introduction

High-quality anti-VASP rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat VASP. Ideal for studying cell migration, adhesion, and platelet activation.

[Learn More](#)

Parameter	Value
Product Name	VASP Rabbit pAb
Remarks/Alias	VASP
Species	Human
GeneID (Human)	7408
GeneID	7408
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human VASP (NP_003361.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50552
Protein Weight	40kDa
Shipping	Ice bag