


Anti-Vav1 Rabbit Polyclonal Antibody For Wb, Elisa - P15498

Item Number: CM0013715



Introduction

This VAV1 rabbit polyclonal antibody is validated for WB and ELISA, and shows cross-reactivity with human, mouse, and rat. Ideal for studying Rho/Rac GTPase signaling in hematopoietic cells and cancer research.

[Learn More](#)

Parameter	Value
Product Name	VAV1 Rabbit pAb
Remarks/Alias	VAV; VAV1
Species	Human
GeneID	7409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 374-590 of human VAV1 (NP_005419.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15498
Protein Weight	98kDa
Shipping	Ice bag