

Anti-Vav2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P52735

Item Number: CM0009880



Introduction

High-specificity rabbit monoclonal antibody targeting human VAV2 (P52735), validated for Western blotting, immunohistochemistry (paraffin), and ELISA. Designed for angiogenesis research, cell migration assays, and vascular biology studies; cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	VAV2 Rabbit mAb
Remarks/Alias	VAV-2; VAV2
Species	Human
GeneID	7410
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human VAV2 (P52735).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52735
Protein Weight	101kDa
Shipping	Ice bag