

Anti-Vav2 Rabbit Polyclonal Antibody For Wb, Elisa - P52735

Item Number: CM0027600



Introduction

VAV2 rabbit polyclonal antibody validated for Western blot and ELISA, recognizing human, mouse, and rat VAV2. Investigate angiogenesis and Rho GTPase signaling with this high-quality antibody.

[Learn More](#)

Parameter	Value
Product Name	VAV2 Rabbit pAb
Remarks/Alias	VAV-2; VAV2
Species	Human
GeneID (Human)	7410
GeneID	7410
Immunogen	Recombinant protein corresponding to amino acids 570-700 of human VAV2 (NP_001127870.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52735
Protein Weight	101kDa
Shipping	Ice bag