


Anti-Vav3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ukw4

Item Number: CM0016604



Introduction

Rabbit polyclonal antibody against VAV3, validated for WB, IF/ICC, and ELISA. Recognizes human and mouse VAV3. Useful for studies on angiogenesis, integrin signaling, and cytoskeleton organization.

[Learn More](#)

Product Name	VAV3 Rabbit pAb
Remarks/Alias	VAV3
Species	Human
GeneID (Human)	10451
GeneID	10451
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 560-680 of human VAV3 (NP_006104.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UKW4
Protein Weight	98kDa
Shipping	Ice bag