


Anti-Vav3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ukw4

Item Number: CM0011925



Introduction

Rabbit monoclonal antibody against human VAV3, validated for WB, IHC-P, ELISA. Detects human and mouse VAV3, a GEF involved in angiogenesis and integrin signaling. Ideal for research in vascular biology and immunology.

[Learn More](#)

Parameter	Value
Product Name	VAV3 Rabbit mAb
Species	Human
Gene ID (Human)	10451
Gene ID	10451
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 700-800 of human VAV3 (Q9UKW4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UKW4
Protein Weight	98kDa
Shipping	Ice bag