


Anti-Vbp1 Polyclonal Antibody For Wb, Elisa - P61758

Item Number: CM0026710



Introduction

High-quality VBP1 Rabbit pAb validated for WB and ELISA. Targets human and mouse Prefoldin subunit 3 (VHL-binding protein 1). Suitable for chaperone and protein folding research.

[Learn More](#)

Parameter	Details
Product Name	VBP1 Rabbit pAb
Remarks/Alias	PFD3; PFDN3; VBP-1; HIBBJ46; VBP1
Species	Human
GeneID (Human)	7411
GeneID	7411
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human VBP1 (NP_003363.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61758
Protein Weight	23kDa
Shipping	Ice bag