

Anti-Vcam-1/Cd106 Rabbit Polyclonal Antibody For Wb, Elisa - P29533

Item Number: CM0019219

Introduction

Rabbit polyclonal antibody against mouse CD106/VCAM-1 for WB and ELISA applications. Specific to mouse, suitable for studies on endothelial cell adhesion, leukocyte migration, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	CD106 Rabbit pAb
Remarks/Alias	CD106; Vcam-1
Species	Mouse
Gene ID	22329
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-310 of mouse CD106. (NP_035823.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P29533
Protein Weight	81kDa
Shipping	Ice bag