

Anti-Vcam1 (Cd106/Incam-100) Polyclonal Antibody For Wb, If-P, Elisa - P19320

Item Number: CM0013821



Introduction

VCAM1 rabbit polyclonal antibody validated for WB, IF-P, ELISA. Reacts with human and mouse VCAM1. Suitable for leukocyte adhesion, inflammation, and immune surveillance research.

[Learn More](#)

Parameter	Value
Product Name	VCAM1 Rabbit pAb
Remarks/Alias	CD106; INCAM-100; VCAM1
Species	Human
GeneID (Human)	7412
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human VCAM1 (NP_001069.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19320
Protein Weight	81kDa
Shipping	Ice bag