


Anti-Vcam1 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P19320 / P29533

Item Number: CM0002985



Introduction

VCAM1 Rabbit PolymAb® validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Supports research on endothelial-leukocyte adhesion and inflammatory signaling pathways. - P19320 / P29533

[Learn More](#)

Parameter	Value
Product Name	VCAM1 Rabbit PolymAb®
Remarks/Alias	CD106; INCAM-100
Species	Human/Mouse
GeneID (Human)	7412
GeneID	7412/22329
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 640-739 of human VCAM1 (P19320).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19320/P29533
Protein Weight	38kDa/71KDa/74KDa/81KDa
Shipping	Ice bag