


Anti-Vcp Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P55072

Item Number: CM0003286



Introduction

High-quality rabbit monoclonal antibody against human VCP (p97) validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for investigating ERAD, Golgi dynamics, DNA repair, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	VCP Rabbit mAb
Remarks/Alias	p97; TERA; CDC48; FTDALS6; VCP
Species	Human
Gene ID (Human)	7415
Gene ID	7415
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VCP (P55072)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55072
Protein Weight	89kDa
Shipping	Ice bag