


Anti-Vcp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P55072

Item Number: CM0018546



Introduction

High-quality rabbit polyclonal antibody targeting human VCP/p97 (TERA) for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse, rat. Ideal for studies on ERAD, autophagy, and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	VCP Rabbit pAb
Remarks/Alias	p97; TERA; CDC48; FTDALS6; VCP
Species	Human
GeneID (Human)	7415
GeneID	7415
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 647-806 of human VCP (NP_009057.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55072
Protein Weight	89kDa
Shipping	Ice bag