

Anti-Vdac 1/2/3 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P21796 / P45880 / Q9Y277

Item Number: CM0005419



Introduction

Recombinant rabbit monoclonal antibody targeting human VDAC1, VDAC2, and VDAC3. Validated for WB, IF-P, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for mitochondrial and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	VDAC 1/2/3 Rabbit PolymAb®
Remarks/Alias	PORIN; VDAC-1
Species	Human
GeneID	7416/7417/7419
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VDAC 1/2/3 (NP_003366.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21796/P45880/Q9Y277
Protein Weight	31kDa
Shipping	Ice bag