

Anti-Vdac1 / Porin Polyclonal Antibody For Wb, Ihc-P, Elisa - P21796

Item Number: CM0020750



Introduction

VDAC1/Porin rabbit polyclonal antibody for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat VDAC1. Ideal for studying mitochondrial ion transport, apoptosis, and cell volume regulation.

[Learn More](#)

Parameter	Value
Product Name	VDAC1 / Porin Rabbit pAb
Remarks/Alias	PORIN; VDAC-1; VDAC1 / Porin
Species	Human
GeneID (Human)	7416
GeneID	7416
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human VDAC1 / Porin (NP_003365.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21796
Protein Weight	31kDa
Shipping	Ice bag