

Anti-Vdac1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P21796

Item Number: CM0003743



Introduction

VDAC1 Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA applications, with human, mouse, rat cross-reactivity. Targets mitochondrial voltage-gated ion channel involved in apoptosis and metabolite transport.

[Learn More](#)

Parameter	Value
Product Name	VDAC1 Rabbit mAb
Remarks/Alias	PORIN; VDAC-1; VDAC1
Species	Human
GeneID (Human)	7416
GeneID	7416
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21796
Protein Weight	31kDa
Shipping	Ice bag