


Anti-Vdac2 Polyclonal Antibody For Wb, Elisa - P45880

Item Number: CM0023640



Introduction

Rabbit polyclonal antibody targeting human VDAC2 (Porin 2), validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for mitochondrial outer membrane ion channel and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	VDAC2 Rabbit pAb
Remarks/Alias	POR; VDAC2
Species	Human
GeneID (Human)	7417
GeneID	7417
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human VDAC2 (NP_001171712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45880
Protein Weight	32kDa
Shipping	Ice bag