

Anti-Vdac2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P45880

Item Number: CM0003492



Introduction

Rabbit monoclonal antibody against human VDAC2, validated for Western blot, IHC-P, and ELISA. Detects endogenous VDAC2 in human, mouse, and rat samples. Suitable for mitochondrial and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	VDAC2 Rabbit mAb
Remarks/Alias	POR; VDAC2
Species	Human
GeneID (Human)	7417
GeneID	7417
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human VDAC2 (NP_003366.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45880
Protein Weight	32kDa
Shipping	Ice bag