


Anti-Vdac2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P45880

Item Number: CM0024368



Introduction

VDAC2 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Recognizes non-selective voltage-gated ion channel VDAC2, a mitochondrial outer membrane porin involved in ion transport and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	VDAC2 Rabbit pAb
Remarks/Alias	POR; VDAC2
Species	Human
GeneID (Human)	7417
GeneID	7417
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VDAC2 (NP_003366.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45880
Protein Weight	32kDa
Shipping	Ice bag