

Anti-Vdac3 (Kd Validated) Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Y277

Item Number: CM0012691

Introduction

VDAC3 rabbit monoclonal antibody (KD validated) for WB, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Targets outer mitochondrial membrane ion channel involved in transport and male fertility. SWISS: Q9Y277.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] VDAC3 Rabbit mAb
Remarks/Alias	VDAC-3; HD-VDAC3
Species	Human
Gene ID (Human)	7419
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human VDAC3 (NP_005653.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y277
Protein Weight	31kDa
Shipping	Ice bag