

Anti-Vdac3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y277

Item Number: CM0015611

Introduction

Rabbit anti-VDAC3 polyclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Detects human and mouse VDAC3; ideal for mitochondrial function and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	VDAC3 Rabbit pAb
Remarks/Alias	VDAC-3; HD-VDAC3; VDAC3
Species	Human
GeneID (Human)	7419
GeneID	7419
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-283 of human VDAC3 (NP_005653.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y277
Protein Weight	31kDa
Shipping	Ice bag