


Anti-Vdr Rabbit Monoclonal Antibody For Wb, Elisa - P11473

Item Number: CM0003946



Introduction

Recombinant rabbit monoclonal antibody targeting human VDR, validated for Western blotting and ELISA. Cross-reacts with mouse. Useful for studying vitamin D3 signaling, calcium homeostasis, and bile acid receptor function.

[Learn More](#)

Parameter	Value
Product Name	VDR Rabbit mAb
Remarks/Alias	NR1I1; PPP1R163; VDR
Species	Human
GeneID (Human)	7421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human VDR (NP_000367.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11473
Protein Weight	48kDa
Shipping	Ice bag