


Anti-Vdr Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Ihc-P, Elisa - P11473

Item Number: CM0004245



Introduction

VDR Rabbit monoclonal antibody validated in WB, IP, IF/ICC, IF-P, IHC-P, ELISA. Detects human, mouse, rat vitamin D3 receptor (~48kDa). Ideal for calcium homeostasis, transcription regulation, and bile acid signaling studies.

[Learn More](#)

Parameter	Value
Product Name	VDR Rabbit mAb
Remarks/Alias	NR1I1; PPP1R163; VDR
Species	Human
GeneID	7421
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11473
Protein Weight	48kDa
Shipping	Ice bag