


Anti-Vdr Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P11473

Item Number: CM0021551



Introduction

High-quality rabbit polyclonal antibody targeting human VDR (Vitamin D3 receptor), validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for studying vitamin D3-mediated transcription and calcium homeostasis.

[Learn More](#)

Parameter	Value
Product Name	VDR Rabbit pAb
Remarks/Alias	NR1I1; PPP1R163; VDR
Species	Human
GeneID	7421
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VDR (NP_000367.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11473
Protein Weight	48kDa
Shipping	Ice bag