

Anti-Ve Cadherin Polyclonal Antibody For Wb, Ihc-P, Elisa - P33151

Item Number: CM0015601

Introduction

High-quality rabbit polyclonal antibody targeting human VE-cadherin (CDH5, CD144) for WB, IHC-P, ELISA applications. Detects human and mouse VE-cadherin, suitable for endothelial cell and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	VE Cadherin Rabbit pAb
Remarks/Alias	7B4; CD144; VE Cadherin
Species	Human
GeneID (Human)	1003
GeneID	1003
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 621-784 of human VE Cadherin (NP_001786.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P33151
Protein Weight	88kDa
Shipping	Ice bag