

Anti-Ve Cadherin Rabbit Polyclonal Antibody For Wb, Elisa - P33151

Item Number: CM0018698

Introduction

VE Cadherin Rabbit pAb is a polyclonal antibody raised in rabbit, validated for WB and ELISA. Reacts with human and mouse VE Cadherin (CD144), a key endothelial adhesion protein (SWISS P33151).

[Learn More](#)

Parameter	Value
Product Name	VE Cadherin Rabbit pAb
Remarks/Alias	7B4; CD144; VE Cadherin
Species	Human
GenelD (Human)	1003
GenelD	1003
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 259-358 of human VE Cadherin (NP_001786.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P33151
Protein Weight	88kDa
Shipping	Ice bag