

Anti-Ve-Cadherin Rabbit Polymab® Monoclonal Antibody For Wb, Ihc-P, Elisa - P33151 / P55284

Item Number: CM0002778

Introduction

High-quality rabbit recombinant monoclonal antibody targeting VE-Cadherin (CDH5). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Supports endothelial cell biology and vascular research.

[Learn More](#)

Parameter	Value
Product Name	VE-Cadherin Rabbit PolymAb®
Remarks/Alias	7B4; CD144
Species	Human
Gene ID (Human)	1003
Gene ID	1003/12562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 259-358 of human VE Cadherin (NP_001786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33151/P55284
Protein Weight	88kDa
Shipping	Ice bag