


Anti-Vegf-D Rabbit Monoclonal Antibody For Wb, Elisa - 043915

Item Number: CM0003564



Introduction

Recombinant rabbit monoclonal antibody targeting human VEGF-D (FIGF) for Western blotting and ELISA. It cross-reacts with mouse and rat, making it suitable for angiogenesis, lymphangiogenesis, and vascular permeability studies across species.

[Learn More](#)

Parameter	Value
Product Name	VEGF-D Rabbit mAb
Remarks/Alias	FIGF; VEGF-D
Species	Human
GeneID (human)	2277
GeneID	2277
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human VEGF-D (043915)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43915
Protein Weight	40kDa
Shipping	Ice bag