


Anti-Vegf-D Rabbit Polyclonal Antibody For Wb And Elisa - O43915

Item Number: CM0017929



Introduction

VEGF-D Rabbit Polyclonal Antibody targets human vascular endothelial growth factor D (FIGF), ideal for WB and ELISA. This antibody supports angiogenesis and lymphangiogenesis research, manufactured for consistent performance.

[Learn More](#)

Parameter	Value
Product Name	VEGF-D Rabbit pAb
Remarks/Alias	FIGF; VEGF-D
Species	Human
GeneID (human)	2277
GeneID	2277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-220 of human VEGF-D (NP_004460.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43915
Protein Weight	40kDa
Shipping	Ice bag